



Percutaneous Burrs

for Forefoot Surgery

Mini-invasive and percutaneous surgery have similar objectives than open surgery thus reducing the size of the incision.

- ▶ In2Bones® supplies a range of 8 percutaneous burrs for mini-invasive and percutaneous surgery.
- ▶ 3 designs.
- ▶ Percutaneous burrs are delivered sterile, by pack of 5, or in an individual packaging, each one in an individual double pouch.

SHANNON STRAIGHT BURRS

Straight Cutting Edge



DIAM 2 - Lg 12 MM



DIAM 2 - Lg 20 MM



DIAM 3.1 - Lg 15 MM



DIAM 4.1 - Lg 15 MM

Helicoidal Cutting Edge



DIAM 2 - Lg 8 MM

Indications:

- Osteotomy of M1 and/or of P1
- All surgical bone procedures of the lesser toes
- Contouring of metatarsal heads

The helicoidal design supplies a cut with extraction of bone chips.

WEDGE CONICAL BURRS



DIAM 2.9 - Lg 13 MM



DIAM 3.9 - Lg 12.75 MM




DIAM 4.3 - Lg 13 MM

Indications:

- Osteotomy of M1 and/or of P1 (Akin)
- Metatarsal osteotomies
- Cheilectomy of Hallux Rigidus (with the 4.3mm burr)

This design supplies a conical cut with extraction of bone chips.



PERCUTANEOUS BURRS

Ref. / SKU	Description	Cutting Length	Diameter	Total length	Packaging	Quantity /Box
G01 00941	Straight burr SHANNON	12 mm	2.0 mm	75 mm	Sterile	5
G01 001511	Straight burr SHANNON	20 mm	2.0 mm	75 mm	Sterile	5
G01 00981	Straight burr SHANNON	15 mm	3.1 mm	75 mm	Sterile	5
G01 00991	Straight burr SHANNON	15 mm	4.1 mm	75 mm	Sterile	5
G01 00951	Straight burr SHANNON	8 mm	2.0 mm	75 mm	Sterile	5
G01 00961	Conical burr WEDGE	13 mm	2.9 mm	75 mm	Sterile	5
G01 00971	Conical burr WEDGE	13 mm	4.3 mm	75 mm	Sterile	5
G01 01451	Conical burr WEDGE	12.75 mm	3.9 mm	75 mm	Sterile	5

INDICATIONS

Burrs are indicated for forefoot mini-invasive and percutaneous surgery.



www.in2bones.com/ifu



RECOMMANDATION

- ▶ It is recommended to carefully read the instructions for use available in the package insert.

DEVICES

- ▶ Instruments connected to a power driver: Classe IIa - CE0086

MANUFACTURER

- ▶ In2Bones SAS
28, chemin du Petit Bois - 69130 Ecully - FRANCE
Phone : +33 (0)4 72 29 26 26 - Fax : +33 (0)4 72 29 26 29

DOCUMENT

- ▶ Reference : BR-DIG-BURRS-EN-062018

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

Always refer to the appropriate instructions for use for complete clinical instructions.

Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.

Except otherwise specified, all commercial trademarks used in this document, In2bones and the In2bones logo are trademarks or registered trademarks of In2bones or its subsidiaries.

